




Copying of this document, and giving it to others and the use or communication of the contents thereof, are forbidden without express authority. Offenders are liable to the payment of damages. All rights are reserved in the event of the grant of a patent or the registration of a utility model or design.

Pratzen- und Stutzenstellung nach Grundriss


Stückliste / Parts list	027654
Benennung / Name	Nummer / Number
Dazugehörige Unterlagen / Additional Documents	



KASAG		KASAG Swiss AG	
Hönggengasse 4, CH-3550 Langnau, Switzerland		TAG Nr.: 01	
Fabrik-Nr.	(y)	Antlagennummer	
Baujahr	2017	Behälter	Mantel
Max. zul. Druck	bar -1 / 6		-1 / 10
Prüfdruck	bar 8,9		16,2
Druckprobe am			
Zul. Min./Max. Temp.	°C -10 / 150		-10 / 150
Volumen	L 1295		60
Medium	Wässrige Lösung		Pyrogenfreies Wasser
Klammerschrauben Typ/Anz.			
Kategorie gemäss 2014/68/EU			IV

Schweissangaben:
 Zusatzwerkstoffe:
 - Produktberührt CN 20/25
 und am Produkt-Mantel / Ki. Böden angeschw. Teile
 - Alle übrigen Thermit GE316L
 - Nicht bezeichnete Schweissnähte 141

Ferritgehalt: produktberührt
 - Grundwerkstoff <3%
 - Schweissnähte <3%






Havariegewicht = 1990 kg	
innen / inside	außen / outside
Welding seams	geschweißte Naht < 0,5 mm
Shells	geschweißte Naht > 0,5 mm
Bolzen / Flanschnuten	Picking / Flanschnuten
Übersichtsskizze	Non
CS Zeichnung	CS Zeichnung


Oberfläche / Surface	
Bezeichnung	Werkstoff
Schweißkonstruktion	Welded construction
Mechanisch bearbeitete Teile	Machined parts
Längengewebe für Halbröhren	Longitudinal measurements for semi-finished material

Ausgangs- und Betriebsdaten	
Design and operating data	
Bezeichnung	Code
Einheit	Unit
Blau- und Druckprüfung	Construction and pressure test
Schweißnahtbehandlung	Preparation of welds
Qualitätsstufe nach EN ISO 9817	Quality level according to EN ISO 9817
Achtung!	Attention!
Bezeichnung	Code
Blau- und Druckprüfung	Construction and pressure test
Schweißnahtbehandlung	Preparation of welds
Qualitätsstufe nach EN ISO 9817	Quality level according to EN ISO 9817
Achtung!	Attention!
Bezeichnung	Code
Einheit	Unit
Blau- und Druckprüfung	Construction and pressure test
Schweißnahtbehandlung	Preparation of welds
Qualitätsstufe nach EN ISO 9817	Quality level according to EN ISO 9817
Achtung!	Attention!

Technische Daten / Technical data							
N01	1	25	ISO	40	EN 1092-1, 11, B1	Test socket insulation	S-E-AT 0149
N02	1	20	ISO	40	EN 1092-1, 11, B1	Inlet tempering media	WN-Flange
N03	1	20	ISO	40	EN 1092-1, 11, B1	Outlet tempering media	WN-Flange
N04	1	B50	-	-	BioControl	Spare (Conductivity)	S-E-AT 0142
N05	1	G1 1/4"	-	-	Ingotl 025 H7	Spare (pH)	S-E-AT 0330
N06	1	40	ISO	40	TC-Connect, ISO 2852	Sampling (9 port)	S-E-AT 0304
N07	1	G 3/8"	-	-	Thermowell TM411	Temperature measurement	S-E-AT 0352
N08	1	25	ISO	40	DIN 11864-3	Product outlet	S-E-AT 0318
N09	1	25	ISO	40	Südmö SVP A386 D-E 30"	Bottom outlet	S-E-AT 0318
N10	1	20	ISO	40	DIN 11864-3-A (BKS)	Inlet pipe (J-tube)	S-E-AT 0349
N11	1	65	ISO	16	DIN 11864-3-A (BKS)	Inlet pipe (Vessel connection)	S-E-AT 0349
N12	1	25	ISO	40	DIN 11864-3-A (BKS)	CIP 2 (Vessel connection)	S-E-AT 0350
N13	1	65	ISO	16	DIN 11864-3-A (BKS)	CIP 2 (Vessel connection)	S-E-AT 0350
N14	1	25	ISO	40	DIN 11864-3-A (BKS)	CIP 1	S-E-AT 0350
N15	1	65	ISO	16	DIN 11864-3-A (BKS)	CIP 1 (Vessel connection)	S-E-AT 0350
N16	1	15	ISO	40	TC-Connect, ISO 2852	Spares port	-
N17	1	40	ISO	40	DIN 11864-3-A (BKS)	Ventilation	-
N18	1	B25	-	-	BioControl	Level switch	S-E-AT 0142
N19	1	40	ISO	16	DIN 32676	Filling level	-
N20	1	B25	-	-	BioControl	Pressure probe	S-E-AT 0143
N21	1	B25	-	-	BioControl	Pressure gauge	S-E-AT 0143
N22	1	32	ISO	40	TC-Connect, ISO 2852	Rupture disc	S-E-AT 0125
N23	1	50	ISO	25	DIN 11864-3-A (BKS)	Sight glass with light	S-E-AT 0348
N24	1	150	16	16	Aseptic block flange	Sight glass	S-E-AT 0106

Stutzen- und Nozzle table							
Stutzen / Nozzle	Material	Material					
N01	1	25	ISO	40	EN 1092-1, 11, B1	Test socket insulation	S-E-AT 0149
N02	1	20	ISO	40	EN 1092-1, 11, B1	Inlet tempering media	WN-Flange
N03	1	20	ISO	40	EN 1092-1, 11, B1	Outlet tempering media	WN-Flange
N04	1	B50	-	-	BioControl	Spare (Conductivity)	S-E-AT 0142
N05	1	G1 1/4"	-	-	Ingotl 025 H7	Spare (pH)	S-E-AT 0330
N06	1	40	ISO	40	TC-Connect, ISO 2852	Sampling (9 port)	S-E-AT 0304
N07	1	G 3/8"	-	-	Thermowell TM411	Temperature measurement	S-E-AT 0352
N08	1	25	ISO	40	DIN 11864-3	Product outlet	S-E-AT 0318
N09	1	25	ISO	40	Südmö SVP A386 D-E 30"	Bottom outlet	S-E-AT 0318
N10	1	20	ISO	40	DIN 11864-3-A (BKS)	Inlet pipe (J-tube)	S-E-AT 0349
N11	1	65	ISO	16	DIN 11864-3-A (BKS)	Inlet pipe (Vessel connection)	S-E-AT 0349
N12	1	25	ISO	40	DIN 11864-3-A (BKS)	CIP 2 (Vessel connection)	S-E-AT 0350
N13	1	65	ISO	16	DIN 11864-3-A (BKS)	CIP 2 (Vessel connection)	S-E-AT 0350
N14	1	25	ISO	40	DIN 11864-3-A (BKS)	CIP 1	S-E-AT 0350
N15	1	65	ISO	16	DIN 11864-3-A (BKS)	CIP 1 (Vessel connection)	S-E-AT 0350
N16	1	15	ISO	40	TC-Connect, ISO 2852	Spares port	-
N17	1	40	ISO	40	DIN 11864-3-A (BKS)	Ventilation	-
N18	1	B25	-	-	BioControl	Level switch	S-E-AT 0142
N19	1	40	ISO	16	DIN 32676	Filling level	-
N20	1	B25	-	-	BioControl	Pressure probe	S-E-AT 0143
N21	1	B25	-	-	BioControl	Pressure gauge	S-E-AT 0143
N22	1	32	ISO	40	TC-Connect, ISO 2852	Rupture disc	S-E-AT 0125
N23	1	50	ISO	25	DIN 11864-3-A (BKS)	Sight glass with light	S-E-AT 0348
N24	1	150	16	16	Aseptic block flange	Sight glass	S-E-AT 0106





Notes:

- 1) Double block and bleed principle applied at the interfaces between CIP and product, media / buffer, clean utilities. Only double block applied for lines which would lose their sterile status by bleeding.
- 2) Number and location of funnels to be defined in next project phase based on 3D design.

Funnel typical bleed lines

2) spool piece to connect a mobile deaerating skid

Master
07.06.17 NST

- 1) 12.06.17 NST
- 2) 11.07.17 THRU
- 3) 13.07.17 NST
- 4) 12.06.17 NST
- 5) 11.07.17 NST
- 6) 11.07.17 NST
- 7) 11.07.17 NST
- 8) 13.07.17 NST
- 9) 13.07.17 NST
- 10) 3.8.17 NST
- 11) 3.8.17 NST
- 12) 3.8.17 NST
- 13) 12.6.17 NST
- 14) 2.8.17 NST
- 15) 3.8.17 NST

X1 12.06.17 NST

TAG-No.	1B4610	1F4611
Name	Buffer Hold Vessel 11	Vent Filter
Design Temperature (°C)	-10-150	-10-150
Design Pressure (bar)	-1/0	-1/0
Material of Construction	1.4539	1.4435
Techn. Data	Volume (l) Vol.-Flow (m³/h)	1000
	Dimensions (mm)	D 1000x H 1750
	Special Features	0.2 µm, 1x5 in
	Producer/Distrib./Year of Const.	
Remarks		

Client: CSL Behring
Project No.: 19004

Architect: ANS
Project No.: ANS Architekten und Planer SIA AG

General contractor: M+W GROUP
Project No.: 230499

Design Partner: M+W GROUP
Project No.: 230499

RELEASED UNDER EPC
ARCHIVING FLAG (DRW CL 2 & 3):
DATE OF IMPRESSION: 201618

CSL Behring AG
RCF Project Lengnau

Revision: 01
Date: 201618
By: [Signature]
Check: [Signature]

Data Sheet Vessel

1B4110 Buffer Hold Vessel 2 (C1 Eq. + Pre/Post Load Wash) 1 (rVIII-SC) / 1UB41 Buffer Hold Chroma 1

Version 06

Status: As Built

This datasheet does also apply to:


2B4750	2B4700	1B4910	1B4610
1B4600	1B4500	1B4300	


Total number: 8

History:


Vers. Date

06.0	22.11.2017
05.0	23.10.2017
04.0	19.12.2016
03.0	08.12.2016
02.0	07.07.2016

Function	Company	Name	Date	Signature
Author	M+W	<i>Aku</i>	22.11.2017	<i>A. Kuhn</i>
Review	<i>M+W</i>	<i>S.M.</i>	<i>21.11.2017</i>	<i>S.M.</i>
Approval				
CSL Behring		M+W		
<p>CSL Behring Biotherapies for Life™</p> <p>CSL Behring Recombinant Facility AG Wankdorfstrasse 10 CH-3000 Bern 22 Switzerland</p>		<p> M+W GROUP</p> <p>M+W Central Europe GmbH Lotterbergstr. 30 D-70499 Stuttgart Germany</p>		
Project Number CSL Behring 16004		Project Number M+W 2304996		
Document Number CSL Behring		Document Number M+W D-P-DA-0077	Version 06.0	
Project RCF Project Lengnau		Document Type / Description Data Sheet	Page 1	

Project-No.		2304996		Data Sheet									
Code		NRCFF		Vessel									
Tag-No.		1B4110											
PFD-No.		PVF_B_01_0060		Building-No.		B		Process		1 (rVIII-SC) / 1UB41 Buffer Hold Chroma 1			
P&ID -No.		PRI_B_01_0071		Level		10		Name		Buffer Hold Vessel 2 (C1 Eq. + Pre/Post Load Wash)			
Drawing-No.		E11047		Room-No.		B_10_2027		Type		Vessel			
01		General								Design Data			
02	5	Inquiry No. / Date		NA /		0		Pressure Vessel Code		AD2000; PED			
03	5	Bid No. / Date		0103705004 / 07.04.2017		3	v	Inside Diameter		1000	mm		
04	5	Order No. / Date		4500971482 / 27.04.2017		5	v	Length w/o Support		1807	mm		
05	5	Standard / Regulation		RS.00042 / RS.00043 / Rs.00044		5	v	Bottom Outlet Height		N/A	mm		
06	5	Inspection		RS.00042 / RS.00043 / Rs.00044		0		Nominal volume		1000	l		
07	5	Manufacturer / Supplier		KASAG / KASAG		5	v	Total volume		1295	l		
08	5	Necessary Certificates		RS.00042 / RS.00043 / Rs.00044				Design Temperature					
09	5	Documentation		RS.00042 / RS.00043 / Rs.00044		4	v	Inside		-10-150	°C		
10	0					4	v	Jacket (Heating / Cooling)		-10-150	°C		
11	0							Design Pressure²					
12		Operating Data						0	v	Inside		-1 / 6	bar
13	3	v	Medium		Process Media		0	v	Jacket (Heating / Cooling)		-1 / 10	bar	
14	3	v	Characteristics		aqueous solution		0	v	Type of bottom		dished end DIN 28011		
15	5	v	Working Volume min./max.		78,2 - 1008	l	5	v	Type of top		removeable, DIN 28011		
16	5	v	Operating Temp. Min./max.		19 - 23	°C			Wall Thickness				
17	5	v	Op. Pressure min./max.²		0 - 2,1	bar	5		Top / Bottom / Cylinder		6 / 6 / 5	mm	
18	5	v	Filling Rate min./max.		N/A	m ³ /h	5		Heating-/ Cooling Jacket		3	mm	
19	5	v	Draining Rate min./max.		N/A	m ³ /h	5		Inliner		N/A	mm	
20	3	v	Density / Bulk Density at [T]		1200	20 kg/m ³ °C	5		Insulation / Insulation Jacket		3	mm	
21	5	v	Specific Heat Capacity			~4.2 kJ/kg K	0		Corrosion Allowance		0	mm	
22	3	v	Dynamic Viscosity at [T]		0.002	20 Pa s °C	5		Welding Factor		acc. PED		
23	0	v	pH-Value min./max.		1 - 13,6		0	v	Vessel Orientation		vertical		
24	5	v	Flash Point		N/A	°C	5		Reinforcing Sheet(s)		no		
25	5	v	Inertisation ²		N/A	mbar	5		Test press. in-/outside²		8.9 / 16.2	bar	
26	0	v	Cleaning in Place		Yes		5		Gaskets / Type		acc. Pipe class		
27	0	v	Medium		0.5M NaOH, 0.1M HNO3		5		Heat Ex. Surface / Content		N/A	m ² / l	
28	0	v	Temperature			<=80 °C			Weight of Vessel				
29	0	v	Sterilisation in Place		Yes		5		Empty / Disaster		720 / 2075	kg	
30	0	v	Medium		pyrogen free steam				Construction Details				
31	0	v	Temperature		<135	°C	0		Heating / Cooling		cylinder		
32	0	v	Heating-/Cooling Medium		Tempering Media		5		Type		coil		
33	0	v	Inlet Temperature		14	°C			Support				
34	0	v	Outlet Temperature		20	°C	5		Type / No. / Norm		brackets / 4 /		
35	0	v	Operating Pressure ²		~3	bar			Fixing				
36	0	v	Density at [T]		1000	25 kg/m ³ °C	5		Type / No. / Norm		lifting lugs / 3 /		
37	0	v	Specific Heat Capacity		4,182	kJ/kg K	5				name plate / 1 /		
38	0	v	Dyn. Viscosity at [T]		0,001	25 Pa s °C	0				Earthing Connector/ 1 /		
39	5	v	Thermal Output (max)		N/A	kW	0				/ /		
40	5	v	Thermal Input (max)		N/A	kW	0		Accessories		/ /		
41	5	v	Heating-/ Cooling Rate		N/A / N/A	°C/min	0		Type / No. / Norm		/ /		
42	0	v	Insulation		yes		0				/ /		
43			Materials						0		/ /		
44	0	v	Product Contacted Parts		1.4539		0	v	Agitator seal				
45	3	v	d-Ferrite Content		Fe <3%		0	v	Arrangement		none		
46	5	v	Gaskets		EPDM peroxid cured		0	v	Aseptic Design		yes		
47	0	v	Sight Glasses		DIN 7080		0						
48	0	v	Inliner		N/A				Surface Treatment				
49	5	v	Non Prod. Contacted Parts / Insulation Jacket		coil:1.4404 / 1.4435 rest:1.4301 / 1.4404		5		Outer surface				
50	5	v	Gaskets		N/A		0		Surface finish		grinded		
51	5	v	Supports		1.4301 / 1.4404		5		Surface Roughness		RA <=1.2µm		
52	5	v	Insulation		Fabr. ISOVER, AGI Q132				Welding Seam		grinded		
53	5	v	Screws, Nuts, Bolts		A4-70; A4		4	v	Inner surface				
54	5	v	Exterior coating				0	v	Surface finish		grinded		
55	5	v	Primer		N/A		0	v	Surface properties		RA <=0.6µm		
56	5	v	Final Coating		N/A		0		Welding Seam		grinded		
57	5	v					0						
58	0	v					0						
59			Remarks										
60			1. Lines marked with "v" contain process information										
61			2. Overpressure. Vacuum is marked with a negative sign.										
62													
63	5												
64	5												
65	5												


Distribution of the original or copies of this document, further use or distribution of its content is not allowed without written permission from M+W. Violators will be obligated to pay fees and penalties.

Project-No.	2304996	Data Sheet				
Code	NRCFF					
Tag-No.	1B4110					
PFD-No.	PVF_B_01_0060	Building-No.	B	Process	1 (rVIII-SC) / 1UB41 Buffer Hold Chroma 1	
P&ID -No.	PRI_B_01_0071	Level	10	Name	Buffer Hold Vessel 2 (C1 Eq. + Pre/Post Load Wash)	
Drawing-No.	E11047	Room-No.	B_10_2027	Type	Vessel	

Rev	Table of Nozzles							
	Ident.	No.	DN	PN	Norm	Flange-/Nozzletype	Sealing Face	Service
5	N03	1	150		similar DIN 28117	aseptic block flange, radial	O-ring	0106 - Sight glass
5	N04	1	50		DIN 11864-3 BKS	Dim. DIN 11866-B	O-ring; Form A	0348 - Sight glass with light
5	N05	1	32		Dim. DIN 11866-B	Na-connect	Flat	0125 - Rupture disc
5	N07	1	B25			Neumo BioControl	O-ring	0142 - Pressure gauge
5	N08	1	B25			Neumo BioControl	O-ring	0142 - Pressure probe
6	N09	1	40		DIN 32676	Dim. DIN 11866-B	Flat	0344 - Filling level
5	N11	1	B25			Neumo BioControl	O-ring	0142 - Level switch
5	N12	1	15		DIN 11864-3 BKS	Dim. DIN 11866-B	O-ring; Form A	0351 - Ventilation
6	N13	1	50		Dim. DIN 11866-B	Na-connect	Flat	0304 - Sampling (Spare port)
5	N14	1	65		DIN 11864-3 BKS	Dim. DIN 11866-B	O-ring; Form A	0350 - CIP 1 (vessel-connection)
5	N14.1	1	25		DIN 11864-3 BKS	Dim. DIN 11866-B	O-ring; Form A	0350 - CIP inlet 1
5	N15	1	65		DIN 11864-3 BKS	Dim. DIN 11866-B	O-ring; Form A	0350 - CIP 2 (vessel-connection)
5	N15.1	1	25		DIN 11864-3 BKS	Dim. DIN 11866-B	O-ring; Form A	0350 - CIP inlet 2
5	N19	1	65		DIN 11864-3 BKS	Dim. DIN 11866-B	O-ring; Form A	0316 - Inlet pipe (J-tube)
5	N19.1	1	20"		DIN 11864-3 BKS	Dim. DIN 11866-B	O-ring; Form A	0316 - inlet pipe
1	N50	1	25		Südmo block flang	Type Südmo SVP	O-ring	0318 - Bottom outlet
5	N52	1	3/8"			Thermowell		0352 - Temperature measurement
5	N55	1	40		Dim. DIN 11866-B	Na-connect; Nova Septum	Flat	0304 - Sampling (5 port)
3	N58	1	G 1 1/4"		Ingold	25H7	O-ring	0330 - Spare (pH)
5	N59	1	B50			Neumo Biocontrol	O-ring	0142 - Spare (conductivity)
5	N64	1	20	40	DIN EN 1092-1 11	welding neck flange	Form B1	Outlet Tempering Media
5	N65	1	20	40	DIN EN 1092-1 11	welding neck flange	Form B1	Inlet Tempering Media
1	N80	1	1/4"		supplier standard	socket with thread		0149 - testsocket insulation


Rev	Remarks Nozzles
0	Nozzle typical number: S-E-AT-XXXX(number in Service column)
5	*DN25 for factor IX/VII vessels
0	
0	
0	
0	

Distribution of the original or copies of this document, or further use or distribution of its content is not allowed without written permission from M +W. Violators will be obligated to pay fees and penalties.

Project-No.	2304996	Data Sheet				
Code	NRCFF					
Tag-No.	1B4110					
		Vessel				
PFD-No.	PVF_B_01_0060	Building-No.	B	Process	1 (rVIII-SC) / 1UB41 Buffer Hold Chroma 1	
P&ID -No.	PRI_B_01_0071	Level	10	Name	Buffer Hold Vessel 2 (C1 Eq. + Pre/Post Load Wash)	
Drawing-No.	E11047	Room-No.	B_10_2027	Type	Vessel	

Sketch


											N 14.1	N 15.1	N 19.1
N01	N03	N04	N05	N07	N08	N09	N11	N12	N13	N14	N15	N19	



Drawing Rev. 01

without written permission from M+W, Violators will be obligated to pay fees and penalties.

Distribution of the original or copies of this document, or further use or distribution of its content is not allowed

Project-No.	2304996	Data Sheet			
Code	NRCFF				
Tag-No.	1B4110				
Vessel					
PFD-No.	PVF_B_01_0060	Building-No.	B	Process	1 (rVIII-SC) / 1UB41 Buffer Hold Chroma 1
P&ID -No.	PRI_B_01_0071	Level	10	Name	Buffer Hold Vessel 2 (C1 Eq. + Pre/Post Load Wash)
Drawing-No.	E11047	Room-No.	B_10_2027	Type	Vessel
Additional Information for Equivalent Equipment					

Tag-No.	Description	Process	PFD-No./PID-No./ Drawing-No.	Building-No./ Level/Room-No.
1B4300	Buffer Hold Vessel 6 (Filter Eq+Displacem.+C2 El.)	1 (rVIII-SC) / 1UB43 Buffer Hold BDI Pooling	PVF_B_01_0063 PRI_B_01_0075 E11048	B 10 B_10_2027
1B4500	Buffer Hold Vessel 7 (C2 Equilibration)	1 (rVIII-SC) / 1UB45 Buffer Hold Chroma 2	PVF_B_01_0064 PRI_B_01_0076 E11050	B 10 B_10_2027
1B4600	Buffer Hold Vessel 10 (C3/C4 Flush + HETP 1)	1 (rVIII-SC) / 1UB46 Buffer Hold Chroma 3	PVF_B_01_0066 PRI_B_01_0077 E11048	B 10 B_10_2027
1B4610	Buffer Hold Vessel 11 (C3/C4 Eq+Post Load Wash)	1 (rVIII-SC) / 1UB46 Buffer Hold Chroma 3	PVF_B_01_0066 PRI_B_01_0078 E11050	B 10 B_10_2027
1B4910	Buffer Hold Vessel 13 (C5 Equilibration + Elution)	1 (rVIII-SC) / 1UB49 Buffer Hold Chroma 5	PVF_B_01_0067 PRI_B_01_0080 E11048	B 10 B_10_2027
2B4700	Buffer Hold Vessel 19 (C3 Eluate 2. Dilution)	2 (rIX-FP) / 2UB47 Buffer Hold Chroma 4	PVF_B_02_0066 PRI_B_02_0087 E11049	B 10 B_10_2020
2B4750	Buffer Hold Vessel 24 (C4 Sanitization)	2 (rIX-FP) / 2UB47 Buffer Hold Chroma 4	PVF_B_02_0068 PRI_B_02_0092 E11049	B 10 B_10_2020

Data Sheet Vessel

1B4110 Buffer Hold Vessel 2 (C1 Eq. + Pre/Post Load Wash) 1 (rVIII-SC) / 1UB41 Buffer Hold Chroma 1

Version 04

Status: RFQ

This datasheet does also apply to:

2B4750

2B4700

Total number:

8

1B4910

1B4610

1B4600


1B4500


1B4300

History:


Vers. Date

04.0	12/19/2016
03.0	12/8/2016
02.0	7/7/2016
01.0	6/22/2016

Function	Company	Name	Date	Signature
Author	M+W		12/19/2016	
Review				
Approval				
CSL Behring		M+W		
 Biotherapies for Life™ CSL Behring Recombinant Facility AG Wankdorfstrasse 10 CH-3000 Bern 22 Switzerland		 M+W Central Europe GmbH Lotterbergstr. 30 D-70499 Stuttgart Germany		
Project Number CSL Behring 16004		Project Number M+W 2304996		
Document Number CSL Behring		Document Number M+W D-P-DA-0077		Version 04.0
Project RCF Project Lengnau		Document Type / Description Data Sheet		Page 1

Project-No.		2304996		Data Sheet								
Code		NRCFF		Vessel								
Tag-No.		1B4110										
PFD-No.		PVF_B_01_0060		Building-No.		B		Process		1 (rVIII-SC) / 1UB41 Buffer Hold Chroma 1		
P&ID -No.		PRI_B_01_0071		Level		10		Name		Buffer Hold Vessel 2 (C1 Eq. + Pre/Post Load Wash)		
Drawing-No.				Room-No.		B_10_2027		Type		Vessel		
01		General						Design Data				
02	0	Inquiry No. / Date	/		0		Pressure Vessel Code	AD2000; PED				
03	0	Bid No. / Date	/		3	v	Inside Diameter	1000	mm			
04	0	Order No. / Date	/		3	v	Length w/o Support	1750	mm			
05	0	Standard / Regulation	compliant to technical specification		0	v	Bottom Outlet Height		mm			
06	0	Inspection	compliant to technical specification		0		Nominal volume	1000	l			
07	0	Manufacturer / Supplier	/		3	v	Total volume	1196	l			
08	0	Necessary Certificates	compliant to technical specification			v	Design Temperature					
09	0	Documentation	compliant to docu requirements		4	v	Inside	-10-150	°C			
10	0				4	v	Jacket (Heating / Cooling)	-10-150	°C			
11	0						Design Pressure²					
12		Operating Data				0	v	Inside	-1	/	6	bar
13	3	Medium	Process Media		0	v	Jacket (Heating / Cooling)	-1	/	10	bar	
14	3	Characteristics	aqueous solution		0	v	Type of bottom	dished end DIN 28011				
15	3	Working Volume min./max.	82.8	1000	l	0	Type of top	dished end DIN 28011				
16	2	Operating Temp. Min./max.			°C		Wall Thickness					
17	2	Op. Pressure min./max. ²			bar	0	Top / Bottom / Cylinder	/	/		mm	
18	0	Filling Rate min./max.			m ³ /h	0	Heating-/ Cooling Jacket				mm	
19	0	Draining Rate min./max.			m ³ /h	0	Inliner				mm	
20	3	Density / Bulk Density at [T]	1200	20	kg/m ³ °C	0	Insulation / Insulation Jacket				mm	
21	2	Specific Heat Capacity			kJ/kg K	0	Corrosion Allowance	0	mm			
22	3	Dynamic Viscosity at [T]	0.002	20	Pa s °C	0	Welding Factor					
23	0	pH-Value min./max.	1 - 13.6			0	Vessel Orientation	vertical				
24	0	Flash Point			°C	0	Reinforcing Sheet(s)					
25	0	Inertisation ²			mbar	0	Test press. in-/outside²				bar	
26	0	Cleaning in Place	Yes			0	Gaskets / Type					
27	0	Medium	0.5M NaOH, 0.1M HNO3			0	Heat Ex. Surface / Content				m ² / l	
28	0	Temperature	<=80		°C		Weight of Vessel					
29	0	Sterilisation in Place	Yes			0	Empty / Disaster	/			kg	
30	0	Medium	pyrogen free steam				Construction Details					
31	0	Temperature	<135		°C	0	Heating / Cooling	cylinder				
32	0	Heating-/Cooling Medium	Tempering Media			4	Type	coil or jacket*3				
33	0	Inlet Temperature	14		°C		Support					
34	0	Outlet Temperature	20		°C	2	Type / No. / Norm	brackets / 4/ acc. Typical				
35	0	Operating Pressure ²	~3		bar		Fixing					
36	0	Density at [T]	1000	25	kg/m ³ °C	0	Type / No. / Norm	lifting lugs / /				
37	0	Specific Heat Capacity	4,182		kJ/kg K	0		name plate / 1/ acc. Typical				
38	0	Dyn. Viscosity at [T]	0,001	25	Pa s °C	0		Earthing Connection/ 1/				
39	0	Thermal Output (max)	cf Spec		kW	0		/ /				
40	0	Thermal Input (max)	cf Spec		kW	0	Accessories	/ /				
41	0	Heating-/ Cooling Rate	/ cf Spec		°C/min	0	Type / No. / Norm	/ /				
42	0	Insulation	yes			0		/ /				
43		Materials				0		/ /				
44	0	Product Contacted Parts	1.4539			0	Agitator seal					
45	3	d-Ferrite Content	Fe <3%			0	Arrangement	none				
46	3	Gaskets	EPDM / MVQ-silicone			0	Aseptic Design	yes				
47	0	Sight Glasses	DIN 7080			0						
48	0	Inliner	N/A				Surface Treatment					
49	2	Non Prod. Contacted Parts / Insulation Jacket	ds/coil:1.4404/1.4435/1.4571 rest:1.4301 or equi.*5)			0	Outer surface					
50	0	Surface finish				0	Surface finish	uniform grinding				
51	0	Gaskets	Gylon			0	Surface Roughness	RA <=1.2µm				
52	0	Supports	V4A			4	Welding Seam	polished eg. Scotch bride				
53	0	Insulation	mineral wool, AS quality			v	Inner surface					
54	0	Screws, Nuts, Bolts	A2-70; A4			4	Surface finish	grinded				
55	0	Exterior coating				0	Surface properties	RA <=0.6µm				
56	0	Primer				0	Welding Seam	grinded				
57	0	Final Coating				0						
58	0											
59		Remarks										
60		1. Lines marked with "v" contain process information										
61		2. Overpressure. Vacuum is marked with a negative sign.										
62												
63	2	4. dynamic loads need to be considered acc to regulations										
64	2	5. Cladding optional in offer: 1.4301 and 1.4404/1.4435										
65	2	6. operation data: see document "scope of order"										


Distribution of the original or copies of this document, further use or distribution of its content is not allowed without written permission from M+W. Violators will be obligated to pay fees and penalties.

Project-No.	2304996	Data Sheet				
Code	NRCFF					
Tag-No.	1B4110					
PFD-No.	PVF_B_01_0060	Building-No.	B	Process	1 (rVIII-SC) / 1UB41 Buffer Hold Chroma 1	
P&ID -No.	PRI_B_01_0071	Level	10	Name	Buffer Hold Vessel 2 (C1 Eq. + Pre/Post Load Wash)	
Drawing-No.		Room-No.	B_1O_2027	Type	Vessel	


Rev	Table of Nozzles							
	Ident.	No.	DN	PN	Norm	Flange-/Nozzletype	Sealing Face	Service
1	N01	1	600			Zimmerlin lid	O-ring	0315 - Manhole with safety switch / indicator
1	N03	1	80		similar DIN 28117	aseptic block flange, radial	O-ring	0106 - Sight glass
1	N04	1	50		similar DIN 28117	aseptic block flange, radial	O-ring	0317 - Sight glass with light
1	N05	1	2"		ASTM	Na-connect, radial	Flat, ISO 2852	0125 - Rupture disc
1	N07	1	B50			Neumo BioControl, radial	O-ring	0143 - Pressure gauge
1	N08	1	B50			Neumo BioControl, radial	O-ring	0143 - Pressure probe
1	N09	1	B50			Neumo BioControl, vertical	O-ring	0142 - filling level
1	N11	1	B50			Neumo BioControl, vertical	O-ring	0142 - Level switch
3	N12	1	15		weld in	Gemü B600	N/A	0337 - Ventilation
1	N13	1	2"		ASTM	Na-connect, radial	Flat, ISO 2852	spare port
1	N14	1	B50			Neumo Biocontrol, radial	O-ring	0303 - CIP 1 (vessel-connection)
3	N14.1	1	25	25	DIN 11864-2 BF	Dim. DIN 11866-B	O-ring; Form A	0303 - CIP inlet 1
1	N15	1	B50			Neumo Biocontrol, radial	O-ring	0303 - CIP 2 (vessel-connection)
3	N15.1	1	25	25	DIN 11864-2 BF	Dim. DIN 11866-B	O-ring; Form A	0303 - CIP inlet 2
1	N19	1	B25			Neumo Biocontrol, radial	O-ring	0316 - Inlet pipe (J-tube)
3	N19.1	1	15	25	DIN 11864-2 BF	Dim. DIN 11866-B	O-ring; Form A	0316 - inlet pipe
1	N50	1	25		Südmo block flang	Type Südmo SVP	O-ring	0318 - Bottom outlet
3	N52	1	N/A			Thermowell		0332 - Temperature measurement
3	N55	1	1 1/2"		ASTM	Na-connect; Nova Septum	Flat, ISO 2852	0304 - Sampling (5 port)
3	N58	1	G 1 1/4"		Ingold	25H7	O-ring	0330 - Spare (pH)
3	N59	1	B50			Neumo Biocontrol, radial	O-ring	0143 - Spare (conductivity)
1	N64	1	25	40	DIN EN 1092-1 11	welding neck flange	Form B1	Outlet Tempering Media
1	N65	1	25	40	DIN EN 1092-1 11	welding neck flange	Form B1	Inlet Tempering Media
1	N80	1	1/4"		supplier standard	socket with thread		0149 - testsocket insulation

Rev	Remarks Nozzles
0	Nozzle typical number: S-E-AT-XXXX(number in Service column)
0	
0	
0	
0	
0	

Distribution of the original or copies of this document, or further use or distribution of its content is not allowed without written permission from M +W. Violators will be obligated to pay fees and penalties.

Project-No.	2304996	Data Sheet				
Code	NRCFF					
Tag-No.	1B4110					
		Vessel				
PFD-No.	PVF_B_01_0060	Building-No.	B	Process	1 (rVIII-SC) / 1UB41 Buffer Hold Chroma 1	
P&ID -No.	PRI_B_01_0071	Level	10	Name	Buffer Hold Vessel 2 (C1 Eq. + Pre/Post Load Wash)	
Drawing-No.		Room-No.	B_10_2027	Type	Vessel	
Sketch						


											N 14.1	N 15.1	N 19.1
N01	N03	N04	N05	N07	N08	N09	N11	N12	N13	N14	N15	N19	



Drawing Rev. 01

without written permission from M +W. Violators will be obligated to pay fees and penalties.

Distribution of the original or copies of this document, or further use or distribution of its content is not allowed

Project-No.	2304996	Data Sheet			
Code	NRCFF				
Tag-No.	1B4110				
Vessel					
PFD-No.	PVF_B_01_0060	Building-No.	B	Process	1 (rVIII-SC) / 1UB41 Buffer Hold Chroma 1
P&ID -No.	PRI_B_01_0071	Level	10	Name	Buffer Hold Vessel 2 (C1 Eq. + Pre/Post Load Wash)
Drawing-No.		Room-No.	B_10_2027	Type	Vessel
Additional Information for Equivalent Equipment					

Tag-No.	Description	Process	PFD-No./PID-No./ Drawing-No.	Building-No./ Level/Room-No.
1B4300	Buffer Hold Vessel 6 (Filter Eq+Displacem.+C2 El.)	1 (rVIII-SC) / 1UB43 Buffer Hold BDI Pooling	PVF_B_01_0063 PRI_B_01_0075	B 10 B_10_2027
1B4500	Buffer Hold Vessel 7 (C2 Equilibration)	1 (rVIII-SC) / 1UB45 Buffer Hold Chroma 2	PVF_B_01_0064 PRI_B_01_0076	B 10 B_10_2027
1B4600	Buffer Hold Vessel 10 (C3/C4 Flush + HETP 1)	1 (rVIII-SC) / 1UB46 Buffer Hold Chroma 3	PVF_B_01_0066 PRI_B_01_0077	B 10 B_10_2027
1B4610	Buffer Hold Vessel 11 (C3/C4 Eq+Post Load Wash)	1 (rVIII-SC) / 1UB46 Buffer Hold Chroma 3	PVF_B_01_0066 PRI_B_01_0078	B 10 B_10_2027
1B4910	Buffer Hold Vessel 13 (C5 Equilibration + Elution)	1 (rVIII-SC) / 1UB49 Buffer Hold Chroma 5	PVF_B_01_0067 PRI_B_01_0080	B 10 B_10_2027
2B4700	Buffer Hold Vessel 19 (C3 Eluate 2. Dilution)	2 (rIX-FP) / 2UB47 Buffer Hold Chroma 4	PVF_B_02_0066 PRI_B_02_0087	B 10 B_10_2020
2B4750	Buffer Hold Vessel 24 (C4 Sanitization)	2 (rIX-FP) / 2UB47 Buffer Hold Chroma 4	PVF_B_02_0068 PRI_B_02_0092	B 10 B_10_2020

Data Sheet Vessel

1B4110 Buffer Hold Vessel 2 (C1 Eq. + Pre/Post Load Wash) 1 (rVIII-SC) / 1UB41 Buffer Hold Chroma 1

Version 06

Status: As Built

This datasheet does also apply to:

2B4750

2B4700

1B4600

Total number:

8

1B4910


1B4300

1B4610

History:


Vers. Date

06.0	22.11.2017
05.0	23.10.2017
04.0	19.12.2016
03.0	08.12.2016
02.0	07.07.2016

Function	Company	Name	Date	Signature
Author	M+W	<i>Aku</i>	22.11.2017	<i>A. Fuhler</i>
Review	<i>M+W</i>	<i>S. M.</i>	<i>22.11.2017</i>	<i>[Signature]</i>
Approval	<i>CSL</i>	<i>Doi</i>	<i>14.12.17</i>	<i>[Signature]</i>
CSL Behring		M+W		
<p>CSL Behring Biotherapies for Life™</p> <p>CSL Behring Recombinant Facility AG Wankdorfstrasse 10 CH-3000 Bern 22 Switzerland</p>		<p> M+W GROUP</p> <p>M+W Central Europe GmbH Lotterbergstr. 30 D-70499 Stuttgart Germany</p>		
Project Number CSL Behring 16004		Project Number M+W 2304996		
Document Number CSL Behring		Document Number M+W D-P-DA-0077		Version 06.0
Project RCF Project Lengnau		Document Type / Description Data Sheet		Page 1

Project-No.		2304996		Data Sheet			M+W GROUP							
Code		NRCFF		Vessel										
Tag-No.		1B4110												
PFD-No.		PVF_B_01_0060		Building-No.		B		Process		1 (rVIII-SC) / 1UB41 Buffer Hold Chroma 1				
P&ID -No.		PRI_B_01_0071		Level		10		Name				Buffer Hold Vessel 2 (C1 Eq. + Pre/Post Load Wash)		
Drawing-No.		E11047		Room-No.		B_10_2027		Type				Vessel		
01	v	General						v	Design Data					
02	5	Inquiry No. / Date		NA /		0	Pressure Vessel Code		AD2000; PED					
03	5	Bid No. / Date		0103705004 / 07.04.2017		3	Inside Diameter		1000 mm					
04	5	Order No. / Date		4500971482 / 27.04.2017		5	Length w/o Support		1807 mm					
05	5	Standard / Regulation		RS.00042 / RS.00043 / Rs.00044		5	Bottom Outlet Height		N/A mm					
06	5	Inspection		RS.00042 / RS.00043 / Rs.00044		0	Nominal volume		1000 l					
07	5	Manufacturer / Supplier		KASAG / KASAG		5	Total volume		1295 l					
08	5	Necessary Certificates		RS.00042 / RS.00043 / Rs.00044		v	Design Temperature							
09	5	Documentation		RS.00042 / RS.00043 / Rs.00044		4	Inside		-10-150		°C			
10	0					4	Jacket (Heating / Cooling)		-10-150		°C			
11	0					v	Design Pressure ²							
12		Operating Data						0	Inside		-1 / 6		bar	
13	3	Medium		Process Media		0	Jacket (Heating / Cooling)		-1 / 10		bar			
14	3	Characteristics		aqueous solution		0	Type of bottom		dished end DIN 28011					
15	5	Working Volume min./max.		78,2 - 1008 l		5	Type of top		removeable, DIN 28011					
16	5	Operating Temp. Min./max.		19 - 23 °C			Wall Thickness							
17	5	Op. Pressure min./max. ²		0 - 2,1 bar		5	Top / Bottom / Cylinder		6 / 6 / 5		mm			
18	5	Filling Rate min./max.		N/A m ³ /h		5	Heating-/ Cooling Jacket		3 mm					
19	5	Draining Rate min./max.		N/A m ³ /h		5	Inliner		N/A mm					
20	3	Density / Bulk Density at [T]		1200 20 kg/m ³ °C		5	Insulation / Insulation Jacket		3 mm					
21	5	Specific Heat Capacity		~4.2 kJ/kg K		0	Corrosion Allowance		0 mm					
22	3	Dynamic Viscosity at [T]		0.002 20 Pa s °C		5	Welding Factor		acc. PED					
23	0	pH-Value min./max.		1 - 13,6		0	Vessel Orientation		vertical					
24	5	Flash Point		N/A °C		5	Reinforcing Sheet(s)		no					
25	5	Inertisation ²		N/A mbar		5	Test press. in-/outside ²		8,9 / 16,2		bar			
26	0	Cleaning in Place		Yes		5	Gaskets / Type		acc. Pipe class					
27	0	Medium		0.5M NaOH, 0.1M HNO3		5	Heat Ex. Surface / Content		N/A		m ² / l			
28	0	Temperature		<=80 °C			Weight of Vessel							
29	0	Sterilisation in Place		Yes		5	Empty / Disaster		720 / 2075		kg			
30	0	Medium		pyrogen free steam			Construction Details							
31	0	Temperature		<135 °C		0	Heating / Cooling		cylinder					
32	0	Heating-/Cooling Medium		Tempering Media		5	Type		coil					
33	0	Inlet Temperature		14 °C			Support							
34	0	Outlet Temperature		20 °C		5	Type / No. / Norm		brackets / 4 /					
35	0	Operating Pressure ²		~3 bar			Fixing							
36	0	Density at [T]		1000 25 kg/m ³ °C		5	Type / No. / Norm		lifting lugs / 3 /					
37	0	Specific Heat Capacity		4,182 kJ/kg K		5			name plate / 1 /					
38	0	Dyn. Viscosity at [T]		0.001 25 Pa s °C		0			earthing Connector/ 1 /					
39	5	Thermal Output (max)		N/A kW		0			/ /					
40	5	Thermal Input (max)		N/A kW		0	Accessories		/ / /					
41	5	Heating-/ Cooling Rate		N/A / N/A °C/min		0	Type / No. / Norm		/ / /					
42	0	Insulation		yes		0			/ / /					
43		Materials						0						
44	0	Product Contacted Parts		1.4539		0	Agitator seal							
45	3	d-Ferrite Content		Fe <3%		0	Arrangement		none					
46	5	Gaskets		EPDM peroxid cured		0	Aseptic Design		yes					
47	0	Sight Glasses		DIN 7080		0								
48	0	Inliner		N/A			Surface Treatment							
49	5	Non Prod. Contacted Parts / Insulation Jacket		coil: 1.4404 / 1.4435 rest: 1.4301 / 1.4404		5	Outer surface		grinded					
50	5	Gaskets		N/A		0	Surface Roughness		RA <=1.2µm					
51	5	Supports		1.4301 / 1.4404		5	Welding Seam		grinded					
52	5	Insulation		Fabr. ISOVER, AGI Q132		v	Inner surface							
53	5	Screws, Nuts, Bolts		A4-70; A4		4	Surface finish		grinded					
54	5	Exterior coating				0	Surface properties		RA <=0.6µm					
55	5	Primer		N/A		0	Welding Seam		grinded					
56	5	Final Coating		N/A		0								
57	5													
58	0													
59		Remarks												
60		1. Lines marked with "v" contain process information												
61		2. Overpressure, Vacuum is marked with a negative sign.												
62														
63	5													
64	5													
65	5													


Distribution of the original or copies of this document, further use or distribution of its content is not allowed without written permission from M+W. Violators will be obligated to pay fees and penalties.

Project-No.	2304996	Data Sheet				
Code	NRCFF					
Tag-No.	1B4110					
Vessel		Building-No.	B	Process	1 (rVIII-SC) / 1UB41 Buffer Hold Chroma 1	
PFD-No.	PVF_B_01_0060	Level	10	Name	Buffer Hold Vessel 2 (C1 Eq. + Pre/Post Load Wash)	
Drawing-No.	E11047	Room-No.	B_10_2027	Type	Vessel	


Rev	Table of Nozzles							
	Ident.	No.	DN	PN	Norm	Flange-/Nozzle type	Sealing Face	Service
5	N03	1	150		similar DIN 28117	aseptic block flange, radial	O-ring	0106 - Sight glass
5	N04	1	50		DIN 11864-3 BKS	Dim. DIN 11866-B	O-ring; Form A	0348 - Sight glass with light
5	N05	1	32		Dim. DIN 11866-B	Na-connect	Flat	0125 - Rupture disc
5	N07	1	B25			Neumo BioControl	O-ring	0142 - Pressure gauge
5	N08	1	B25			Neumo BioControl	O-ring	0142 - Pressure probe
6	N09	1	40		DIN 32676	Dim. DIN 11866-B	Flat	0344 - Filling level
5	N11	1	B25			Neumo BioControl	O-ring	0142 - Level switch
5	N12	1	15		DIN 11864-3 BKS	Dim. DIN 11866-B	O-ring; Form A	0351 - Ventilation
6	N13	1	50		Dim. DIN 11866-B	Na-connect	Flat	0304 - Sampling (Spare port)
5	N14	1	65		DIN 11864-3 BKS	Dim. DIN 11866-B	O-ring; Form A	0350 - CIP 1 (vessel-connection)
5	N14.1	1	25		DIN 11864-3 BKS	Dim. DIN 11866-B	O-ring; Form A	0350 - CIP inlet 1
5	N15	1	65		DIN 11864-3 BKS	Dim. DIN 11866-B	O-ring; Form A	0350 - CIP 2 (vessel-connection)
5	N15.1	1	25		DIN 11864-3 BKS	Dim. DIN 11866-B	O-ring; Form A	0350 - CIP inlet 2
5	N19	1	65		DIN 11864-3 BKS	Dim. DIN 11866-B	O-ring; Form A	0316 - Inlet pipe (J-tube)
5	N19.1	1	20"		DIN 11864-3 BKS	Dim. DIN 11866-B	O-ring; Form A	0316 - inlet pipe
1	N50	1	25		Südmo block flang	Type Südmo SVP	O-ring	0318 - Bottom outlet
5	N52	1	3/8"			Thermowell		0352 - Temperature measurement
5	N55	1	40		Dim. DIN 11866-B	Na-connect; Nova Septum	Flat	0304 - Sampling (5 port)
3	N58	1	G 1 1/4"		Ingold	25H7	O-ring	0330 - Spare (pH)
5	N59	1	B50			Neumo Biocontrol	O-ring	0142 - Spare (conductivity)
5	N64	1	20	40	DIN EN 1092-1 11	welding neck flange	Form B1	Outlet Tempering Media
5	N65	1	20	40	DIN EN 1092-1 11	welding neck flange	Form B1	Inlet Tempering Media
1	N80	1	1/4"		supplier standard	socket with thread		0149 - testsocket insulation


Distribution of the original or copies of this document, or further use or distribution of its content is not allowed without written permission from M+W. Violators will be obligated to pay fees and penalties.

Rev	Remarks Nozzles
0	Nozzle typical number: S-E-AT-XXXX(number in Service column)
5	*DN25 for factor IX/VII vessels
0	
0	
0	
0	
0	

Project-No.	2304996	Data Sheet			
Code	NRCFF				
Tag-No.	1B4110				
Vessel		Building-No.	B	Process	1 (rVIII-SC) / 1UB41 Buffer Hold Chroma 1
P&ID -No.	PVF_B_01_0060	Level	10	Name	Buffer Hold Vessel 2 (C1 Eq. + Pre/Post Load Wash)
Drawing-No.	E11047	Room-No.	B_10_2027	Type	Vessel
Sketch					

Distribution of the original or copies of this document, or further use or distribution of its content is not allowed without written permission from M+W. Violators will be obligated to pay fees and penalties.



Project-No.	2304996	Data Sheet			
Code	NRCFF				
Tag-No.	1B4110	Vessel			
PFD-No.	PVF_B_01_0060	Building-No.	B	Process	1 (rVIII-SC) / 1UB41 Buffer Hold Chroma 1
P&ID -No.	PRI_B_01_0071	Level	10	Name	Buffer Hold Vessel 2 (C1 Eq. + Pre/Post Load Wash)
Drawing-No.	E11047	Room-No.	B_10_2027	Type	Vessel
Additional Information for Equivalent Equipment					

Tag-No.	Description	Process	PFD-No./PID-No./ Drawing-No.	Building-No./ Level/Room-No.
1B4300	Buffer Hold Vessel 6 (Filter Eq+Displacem.+C2 El.)	1 (rVIII-SC) / 1UB43 Buffer Hold BDI Pooling	PVF_B_01_0063 PRI_B_01_0075 E11048	B 10 B_10_2027
1B4500	Buffer Hold Vessel 7 (C2 Equilibration)	1 (rVIII-SC) / 1UB45 Buffer Hold Chroma 2	PVF_B_01_0064 PRI_B_01_0076 E11050	B 10 B_10_2027
1B4600	Buffer Hold Vessel 10 (C3/C4 Flush + HETP 1)	1 (rVIII-SC) / 1UB46 Buffer Hold Chroma 3	PVF_B_01_0066 PRI_B_01_0077 E11048	B 10 B_10_2027
1B4610	Buffer Hold Vessel 11 (C3/C4 Eq+Post Load Wash)	1 (rVIII-SC) / 1UB46 Buffer Hold Chroma 3	PVF_B_01_0066 PRI_B_01_0078 E11050	B 10 B_10_2027
1B4910	Buffer Hold Vessel 13 (C5 Equilibration + Elution)	1 (rVIII-SC) / 1UB49 Buffer Hold Chroma 5	PVF_B_01_0067 PRI_B_01_0080 E11048	B 10 B_10_2027
2B4700	Buffer Hold Vessel 19 (C3 Eluate 2. Dilution)	2 (rIX-FP) / 2UB47 Buffer Hold Chroma 4	PVF_B_02_0066 PRI_B_02_0087 E11049	B 10 B_10_2020
2B4750	Buffer Hold Vessel 24 (C4 Sanitization)	2 (rIX-FP) / 2UB47 Buffer Hold Chroma 4	PVF_B_02_0068 PRI_B_02_0092 E11049	B 10 B_10_2020

